

Laminin subunit gamma-1

Catalog No: #AP84298



Package Size: #AP84298-1 50ug #AP84298-2 100ug #AP84298-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Laminin subunit gamma-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Laminin B2 chain,Laminin-1 subunit gamma,Laminin-10 subunit gamma,Laminin-11 subunit gamma,Laminin-2 subunit gamma,Laminin-3 subunit gamma,Laminin-4 subunit gamma,Laminin-6 subunit gamma,Laminin-7 subunit gamma,Laminin-8 subunit gamma,Laminin-9 subunit gamma,S-laminin subunit gamma,S-LAM gamma,Lamb-2,Lamc-1
Accession No.	Uniprot ID:P02468
Uniprot	P02468
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lamc1

Note: This product is for in vitro research use only