

## Hepatocyte nuclear factor 1-alpha

Catalog No: #AP84299



Package Size: #AP84299-1 50ug #AP84299-2 100ug #AP84299-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Hepatocyte nuclear factor 1-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	HNF-1-alpha,Liver-specific transcription factor LF-B1,LFB1,Transcription factor 1,TCF-1,Hnf-1,Hnf-1a,Tcf1
Accession No.	Uniprot ID:P15257
Uniprot	P15257
GeneID	24817;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Hnf1a

Note: This product is for in vitro research use only