

Heat shock 70 kDa protein 12A

Catalog No: #AP84302



Package Size: #AP84302-1 50ug #AP84302-2 100ug #AP84302-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Heat shock 70 kDa protein 12A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Kiaa0417
Accession No.	Uniprot ID:Q8K0U4
Uniprot	Q8K0U4
GeneID	73442;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Hspa12a

Note: This product is for in vitro research use only