

Interleukin-34

Catalog No: #AP84304



Package Size: #AP84304-1 50ug #AP84304-2 100ug #AP84304-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-34
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	IL-34
Accession No.	Uniprot ID:Q4KM46
Uniprot	Q4KM46
GeneID	498951;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IL34

Note: This product is for in vitro research use only