

Interleukin-36 gamma

Catalog No: #AP84305



Package Size: #AP84305-1 50ug #AP84305-2 100ug #AP84305-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-36 gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-1-related protein 2,IL-1RP2,Interleukin-1 epsilon,IL-1 epsilon,Interleukin-1 family member 9,IL-1F9,Interleukin-1 homolog 1,IL-1H1,IL1E,IL1F9,IL1H1,IL1RP2,UNQ2456/PRO5737
Accession No.	Uniprot ID:Q9NZH8
Uniprot	Q9NZH8
GeneID	56300;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL36G

Note: This product is for in vitro research use only