

Insulin-like peptide INSL5

Catalog No: #AP84306



Package Size: #AP84306-1 50ug #AP84306-2 100ug #AP84306-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin-like peptide INSL5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Insulin-like peptide 5,Relaxin/insulin-like factor 2,Relaxin/insulin-like protein,Rif,Rif2,Zins3
Accession No.	Uniprot ID:Q9WUG6
Uniprot	Q9WUG6
GeneID	23919;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Insl5

Note: This product is for in vitro research use only