

## Type I iodothyronine deiodinase

Catalog No: #AP84307



Package Size: #AP84307-1 50ug #AP84307-2 100ug #AP84307-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Type I iodothyronine deiodinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	5DI,DIOI,Type 1 DI,Type-I 5'-deiodinase,Itdi1,Txd1
Accession No.	Uniprot ID:Q61153
Uniprot	Q61153
GeneID	13370;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Dio1

Note: This product is for in vitro research use only