

Iron-responsive element-binding protein 2

Catalog No: #AP84309



Package Size: #AP84309-1 50ug #AP84309-2 100ug #AP84309-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Iron-responsive element-binding protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	IRE-BP 2,RCJMB04_5d21
Accession No.	Uniprot ID:Q5ZLQ4
Uniprot	Q5ZLQ4
GeneID	427490;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IREB2

Note: This product is for in vitro research use only