

Insulin receptor substrate 4

Catalog No: #AP84310



Package Size: #AP84310-1 50ug #AP84310-2 100ug #AP84310-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin receptor substrate 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IRS-4, Phosphoprotein of 160 kDa, pp160
Accession No.	Uniprot ID:Q9Z0Y7
Uniprot	Q9Z0Y7
GeneID	16370;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Irs4

Note: This product is for in vitro research use only