

ATP-sensitive inward rectifier potassium channel 10

Catalog No: #AP84312



Package Size: #AP84312-1 50ug #AP84312-2 100ug #AP84312-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ATP-sensitive inward rectifier potassium channel 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ATP-dependent inwardly rectifying potassium channel Kir4.1, Inward rectifier K(+) channel Kir1.2, Potassium channel, inwardly rectifying subfamily J member 10
Accession No.	Uniprot ID:P78508
Uniprot	P78508
GeneID	3766;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: KCNJ10

Note: This product is for in vitro research use only