

Integrin alpha-5

Catalog No: #AP84313



Package Size: #AP84313-1 50ug #AP84313-2 100ug #AP84313-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Integrin alpha-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CD49 antigen-like family member E,Fibronectin receptor subunit alpha,Integrin alpha-F,VLA-5,CD49e
Accession No.	Uniprot ID:P11688
Uniprot	P11688
GeneID	16402;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Itga5

Note: This product is for in vitro research use only