

Inter-alpha-trypsin inhibitor heavy chain H6

Catalog No: #AP84314



Package Size: #AP84314-1 50ug #AP84314-2 100ug #AP84314-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Inter-alpha-trypsin inhibitor heavy chain H6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Inter-alpha-trypsin inhibitor heavy chain H5-like protein, Inter-alpha inhibitor H5-like protein, ITIH5L, UNQ6369/PRO21074
Accession No.	Uniprot ID:Q6UXX5
Uniprot	Q6UXX5
GeneID	347365;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ITIH6

Note: This product is for in vitro research use only