

Heterogeneous nuclear ribonucleoprotein F

Catalog No: #AP84315



Package Size: #AP84315-1 50ug #AP84315-2 100ug #AP84315-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heterogeneous nuclear ribonucleoprotein F
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	hnRNP F, HNRPF
Accession No.	Uniprot ID:Q5E9J1
Uniprot	Q5E9J1
GeneID	506917;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HNRNPF

Note: This product is for in vitro research use only