

Heat shock protein beta-6

Catalog No: #AP84316



Package Size: #AP84316-1 50ug #AP84316-2 100ug #AP84316-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Heat shock protein beta-6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HspB6,Gm479
Accession No.	Uniprot ID:Q5EBG6
Uniprot	Q5EBG6
GeneID	243912;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hspb6

Note: This product is for in vitro research use only