

Hyaluronidase-1

Catalog No: #AP84317



Package Size: #AP84317-1 50ug #AP84317-2 100ug #AP84317-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hyaluronidase-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Hyal-1,Hyaluronoglucosaminidase-1
Accession No.	Uniprot ID:Q91ZJ9
Uniprot	Q91ZJ9
GeneID	15586;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hyal1

Note: This product is for in vitro research use only