

Interferon alpha-1/13

Catalog No: #AP84318



Package Size: #AP84318-1 50ug #AP84318-2 100ug #AP84318-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interferon alpha-1/13
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IFN-alpha-1/13,Interferon alpha-D,LeIF D
Accession No.	Uniprot ID:P01562
Uniprot	P01562
GeneID	3439;3447;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IFNA1

Note: This product is for in vitro research use only