

Interferon alpha-6

Catalog No: #AP84319



Package Size: #AP84319-1 50ug #AP84319-2 100ug #AP84319-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interferon alpha-6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IFN-alpha-6, Interferon alpha-8, IFN-alpha-8, lfa6, lfa8, lfna8
Accession No.	Uniprot ID:P07350
Uniprot	P07350
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: lfna6

Note: This product is for in vitro research use only