

Iodotyrosine deiodinase 1

Catalog No: #AP84322



Package Size: #AP84322-1 50ug #AP84322-2 100ug #AP84322-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Iodotyrosine deiodinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IYD-1,Iodotyrosine dehalogenase 1,Dehal1
Accession No.	Uniprot ID:Q9DCX8
Uniprot	Q9DCX8
GeneID	70337;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:lyd

Note: This product is for in vitro research use only