

Homocysteine-responsive endoplasmic reticulum-resident ubiquitin-like domain member 1 protein

Catalog No: #AP84323

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #AP84323-1 50ug #AP84323-2 100ug #AP84323-3 1mg

Description

Product Name	Homocysteine-responsive endoplasmic reticulum-resident ubiquitin-like domain member 1 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q9JJK5
Uniprot	Q9JJK5
GeneID	64209;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Herpud1

Note: This product is for in vitro research use only