

Hemoglobin subunit beta-1

Catalog No: #AP84327



Package Size: #AP84327-1 50ug #AP84327-2 100ug #AP84327-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hemoglobin subunit beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Beta-1-globin,Hemoglobin beta chain, major-form,Hemoglobin beta-1 chain
Accession No.	Uniprot ID:P02091
Uniprot	P02091
GeneID	24440;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hbb

Note: This product is for in vitro research use only