

## Hemoglobin subunit gamma-1

Catalog No: #AP84328



Package Size: #AP84328-1 50ug #AP84328-2 100ug #AP84328-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Hemoglobin subunit gamma-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Gamma-1-globin,Hb F Agamma,Hemoglobin gamma-1 chain,Hemoglobin gamma-A chain,PRO2979
Accession No.	Uniprot ID:P69891
Uniprot	P69891
GeneID	3047;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HBG1

Note: This product is for in vitro research use only