

Hypoxia-inducible factor 1-alpha

Catalog No: #AP84330



Package Size: #AP84330-1 50ug #AP84330-2 100ug #AP84330-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hypoxia-inducible factor 1-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HIF-1-alpha,ARNT-interacting protein
Accession No.	Uniprot ID:Q61221
Uniprot	Q61221
GeneID	15251;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hif1a

Note: This product is for in vitro research use only