

Interleukin-3

Catalog No: #AP84331



Package Size: #AP84331-1 50ug #AP84331-2 100ug #AP84331-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-3,Hematopoietic growth factor,Mast cell growth factor,MCGF,Multipotential colony-stimulating factor,P-cell-stimulating factor,Csfmu,IL-3
Accession No.	Uniprot ID:P01586
Uniprot	P01586
GeneID	16187;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL3

Note: This product is for in vitro research use only