

## Insulin receptor substrate 1

Catalog No: #AP84336



Package Size: #AP84336-1 50ug #AP84336-2 100ug #AP84336-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Insulin receptor substrate 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	IRS-1,pp185,Irs-1
Accession No.	Uniprot ID:P35570
Uniprot	P35570
GeneID	25467;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Irs1

Note: This product is for in vitro research use only