

Interleukin-5

Catalog No: #AP84343



Package Size: #AP84343-1 50ug #AP84343-2 100ug #AP84343-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-5,B-cell growth factor II,BCGF-II,Cytotoxic T-lymphocyte inducer,Eosinophil differentiation factor,T-cell replacing factor,TRF,II-5
Accession No.	Uniprot ID:P04401
Uniprot	P04401
GeneID	16191;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL5

Note: This product is for in vitro research use only