

## H-2 class II histocompatibility antigen gamma chain

Catalog No: #AP84347



Package Size: #AP84347-1 50ug #AP84347-2 100ug #AP84347-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	H-2 class II histocompatibility antigen gamma chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Ia antigen-associated invariant chain, Ii, MHC class II-associated invariant chain, CD74
Accession No.	Uniprot ID:P10247
Uniprot	P10247
GeneID	25599;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Cd74

Note: This product is for in vitro research use only