

## Hypoxia-inducible factor 1-alpha inhibitor

Catalog No: #AP84348



Package Size: #AP84348-1 50ug #AP84348-2 100ug #AP84348-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Hypoxia-inducible factor 1-alpha inhibitor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Factor inhibiting HIF-1,FIH-1,Hypoxia-inducible factor asparagine hydroxylase,FIH1
Accession No.	Uniprot ID:Q9NWT6
Uniprot	Q9NWT6
GeneID	55662;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HIF1AN

Note: This product is for in vitro research use only