

High mobility group protein B3

Catalog No: #AP84349



Package Size: #AP84349-1 50ug #AP84349-2 100ug #AP84349-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	High mobility group protein B3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	High mobility group protein 2a,HMG-2a,High mobility group protein 4,HMG-4,HMG2A,HMG4
Accession No.	Uniprot ID:O15347
Uniprot	O15347
GeneID	3149;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HMGB3

Note: This product is for in vitro research use only