

# Interleukin-1 receptor type 1

Catalog No: #AP84351

Package Size: #AP84351-1 50ug #AP84351-2 100ug #AP84351-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-1 receptor type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-1R-1,CD121 antigen-like family member A,Interleukin-1 receptor alpha,IL-1R-alpha,Interleukin-1 receptor type I,p80,CD121a,IL1R,IL1RA,IL1RT1
Accession No.	Uniprot ID:P14778
Uniprot	P14778
GeneID	3554;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IL1R1

Note: This product is for in vitro research use only