

Glucosidase 2 subunit beta

Catalog No: #AP84353



Package Size: #AP84353-1 50ug #AP84353-2 100ug #AP84353-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glucosidase 2 subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	80K-H protein,Glucosidase II subunit beta,Protein kinase C substrate 60.1 kDa protein heavy chain,PKC ζ ,G19P1
Accession No.	Uniprot ID:P14314
Uniprot	P14314
GeneID	5589;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PRKCSH

Note: This product is for in vitro research use only