

Interleukin-23 subunit alpha

Catalog No: #AP84354



Package Size: #AP84354-1 50ug #AP84354-2 100ug #AP84354-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-23 subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	IL-23 subunit alpha, Interleukin-23 subunit p19, IL-23p19, SGRF
Accession No.	Uniprot ID:Q9N2H9
Uniprot	Q9N2H9
GeneID	100155546;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IL23A

Note: This product is for in vitro research use only