

Interleukin-28A

Catalog No: #AP84355



Package Size: #AP84355-1 50ug #AP84355-2 100ug #AP84355-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-28A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-28A,Cytokine Zcyto20,Interferon lambda-2,IFN-lambda-2,IFNL2,ZCYTO20
Accession No.	Uniprot ID:Q8IZJ0
Uniprot	Q8IZJ0
GeneID	282616;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL28A

Note: This product is for in vitro research use only