

## Interleukin-6 receptor subunit beta

Catalog No: #AP84356



Package Size: #AP84356-1 50ug #AP84356-2 100ug #AP84356-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-6 receptor subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-6 receptor subunit beta,CDw130,Interleukin-6 signal transducer,Membrane glycoprotein 130,gp130,Oncostatin-M receptor subunit alpha,CD130
Accession No.	Uniprot ID:P40189
Uniprot	P40189
GeneID	3572;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IL6ST

Note: This product is for in vitro research use only