

NFAT2 (Phospho-Ser172) Antibody

Catalog No: #12620

Package Size: #12620-1 50ul #12620-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	NFAT2 (Phospho-Ser172) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms
Specificity	NFAT2 (Phospho-Ser172) Antibody detects endogenous levels of NFAT2 only when phosphorylated at Ser172
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human NFAT2 (Phospho-Ser172)
Target Name	NFAT2
Modification	Phospho
Other Names	NFATC1, NFATc, NFAT2, NF-ATC, NF-ATC transcription factor, NF-ATc1
Accession No.	Swiss-Prot#: O95644NCBI Gene ID: 4772
Uniprot	O95644
GeneID	4772;
Target Species	human
Calculated MW	78kd,101kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images

Western Blot analysis of lysates of 293T cell, using primary antibody at 1:1000 dilution.



Note: This product is for in vitro research use only