

Heat shock 70 kDa protein 4

Catalog No: #AP84360



Package Size: #AP84360-1 50ug #AP84360-2 100ug #AP84360-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Heat shock 70 kDa protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Heat shock 70-related protein APG-2,App2,Hsp110
Accession No.	Uniprot ID:Q61316
Uniprot	Q61316
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hspa4

Note: This product is for in vitro research use only