

## Inorganic pyrophosphatase

Catalog No: #AP84363



Package Size: #AP84363-1 50ug #AP84363-2 100ug #AP84363-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Inorganic pyrophosphatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Pyrophosphate phospho-hydrolase,PPase,PP
Accession No.	Uniprot ID:P37980
Uniprot	P37980
GeneID	280701;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PPA1

Note: This product is for in vitro research use only