

## Iron-responsive element-binding protein 2

Catalog No: #AP84364



Package Size: #AP84364-1 50ug #AP84364-2 100ug #AP84364-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Iron-responsive element-binding protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IRE-BP 2,Iron regulatory protein 2,IRP2,Irp2
Accession No.	Uniprot ID:Q811J3
Uniprot	Q811J3
GeneID	64602;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ireb2

Note: This product is for in vitro research use only