

Histone deacetylase 8

Catalog No: #AP84365



Package Size: #AP84365-1 50ug #AP84365-2 100ug #AP84365-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histone deacetylase 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HD8,HDACL1,CDA07
Accession No.	Uniprot ID:Q9BY41
Uniprot	Q9BY41
GeneID	55869;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HDAC8

Note: This product is for in vitro research use only