

Histone deacetylase 9

Catalog No: #AP84366



Package Size: #AP84366-1 50ug #AP84366-2 100ug #AP84366-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histone deacetylase 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HD9,Histone deacetylase 7B,HD7,Histone deacetylase-related protein,MEF2-interacting transcription repressor MITR,HDAC7,HDAC7B,HDRP,KIAA0744,MITR
Accession No.	Uniprot ID:Q9UKV0
Uniprot	Q9UKV0
GeneID	9734;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HDAC9

Note: This product is for in vitro research use only