

Interleukin-1 receptor-associated kinase 3

Catalog No: #AP84367



Package Size: #AP84367-1 50ug #AP84367-2 100ug #AP84367-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-1 receptor-associated kinase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IRAK-3,IL-1 receptor-associated kinase M,IRAK-M
Accession No.	Uniprot ID:Q8K4B2
Uniprot	Q8K4B2
GeneID	73914;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Irak3

Note: This product is for in vitro research use only