

Interleukin-3 receptor subunit alpha

Catalog No: #AP84373



Package Size: #AP84373-1 50ug #AP84373-2 100ug #AP84373-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-3 receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-3 receptor subunit alpha,CD123,IL3R
Accession No.	Uniprot ID:P26951
Uniprot	P26951
GeneID	3563;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL3RA

Note: This product is for in vitro research use only