

Interferon gamma receptor 1

Catalog No: #AP84374



Package Size: #AP84374-1 50ug #AP84374-2 100ug #AP84374-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interferon gamma receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IFN-gamma receptor 1,CDw119,Interferon gamma receptor alpha-chain,IFN-gamma-R-alpha,CD119
Accession No.	Uniprot ID:P15260
Uniprot	P15260
GeneID	3459;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IFNGR1

Note: This product is for in vitro research use only