

Holliday junction recognition protein

Catalog No: #AP84375



Package Size: #AP84375-1 50ug #AP84375-2 100ug #AP84375-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Holliday junction recognition protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	14-3-3-associated AKT substrate,Fetal liver-expressing gene 1 protein,Up-regulated in lung cancer 9,FAKTS,FLEG1,URLC9
Accession No.	Uniprot ID:Q8NCD3
Uniprot	Q8NCD3
GeneID	55355;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HJURP

Note: This product is for in vitro research use only