

Heterogeneous nuclear ribonucleoprotein U

Catalog No: #AP84376



Package Size: #AP84376-1 50ug #AP84376-2 100ug #AP84376-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heterogeneous nuclear ribonucleoprotein U
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	hnRNP U,Scaffold attachment factor A,SAF-A,p120,pp120,HNRPU,SAFA,U21.1
Accession No.	Uniprot ID:Q00839
Uniprot	Q00839
GeneID	3192;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HNRNPU

Note: This product is for in vitro research use only