

Interferon omega-1

Catalog No: #AP84377



Package Size: #AP84377-1 50ug #AP84377-2 100ug #AP84377-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interferon omega-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	IFN-omega-c1, Interferon alpha-II-1, IFNW
Accession No.	Uniprot ID:P07352
Uniprot	P07352
GeneID	281847;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IFNW1

Note: This product is for in vitro research use only