

Interleukin-27 subunit alpha

Catalog No: #AP84378



Package Size: #AP84378-1 50ug #AP84378-2 100ug #AP84378-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-27 subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	IL-27 subunit alpha,p28,IL27A
Accession No.	Uniprot ID:Q5S1V9
Uniprot	Q5S1V9
GeneID	493187;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL27

Note: This product is for in vitro research use only