

Insulin-like 3

Catalog No: #AP84379



Package Size: #AP84379-1 50ug #AP84379-2 100ug #AP84379-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin-like 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Leydig insulin-like peptide,Ley-I-L,Relaxin-like factor,Rlf
Accession No.	Uniprot ID:Q9WUK0
Uniprot	Q9WUK0
GeneID	114215;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:InsI3

Note: This product is for in vitro research use only