

Phosphoribosyl pyrophosphate synthase-associated protein 1

Catalog No: #AP84382

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP84382-1 50ug #AP84382-2 100ug #AP84382-3 1mg

Description

Product Name	Phosphoribosyl pyrophosphate synthase-associated protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PRPP synthase-associated protein 1,39 kDa phosphoribosypyrophosphate synthase-associated protein,PAP39
Accession No.	Uniprot ID:Q9D0M1
Uniprot	Q9D0M1
GeneID	67763;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Prpsap1

Note: This product is for in vitro research use only