

6-phosphofructokinase, muscle type

Catalog No: #AP84391



Package Size: #AP84391-1 50ug #AP84391-2 100ug #AP84391-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	6-phosphofructokinase, muscle type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Phosphofructo-1-kinase isozyme A, PFK-A, Phosphofructokinase-M, Phosphofructokinase 1, Phosphohexokinase
Accession No.	Uniprot ID:Q0IIG5
Uniprot	Q0IIG5
GeneID	506544;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PFKM

Note: This product is for in vitro research use only