

Calcium/calmodulin-dependent protein kinase type II subunit alpha



Catalog No: #AP84392

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP84392-1 50ug #AP84392-2 100ug #AP84392-3 1mg

Description

Product Name	Calcium/calmodulin-dependent protein kinase type II subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CaM kinase II subunit alpha
Accession No.	Uniprot ID:P11798
Uniprot	P11798
GeneID	12322;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Camk2a

Note: This product is for in vitro research use only